Bottle Top Dispenser

Features & Benefits Easy Maintenance

The internal glass barrel can be easily removed, allowing quick and thorough cleaning. In the event of a breakage, this allows the barrel to be replaced, saving on long term costs.

Confident Dispensing

• High quality construction materials (borosilicate glass, ETFE, FEP, PVA, PP) offer a comprehensive chemical resistance profile.

Adapters

• Supplied with a range of adapters to suit most common laboratory bottles.

Autoclavable

• Fully autoclavable, ensuring sterility.

Range


• Four different volume ranges from 0.5mL up to 50mL.



qe

Liquid Handling

Bottle Top Dispenser - Operation



Liquid Handling

Bottle Top Dispenser

Please note: information provided in the below Chemical Resistance Guide is a guide only and does not represent a guarantee. The information refers to the following construction materials only; borosilicate glass, PTFE, FEP.



Resistant Partially Resistant Non-Resistant

Chemical Resistance Guide				
Chemical Name	Resistance	Chemical Name Resistance		
1,4 Dioxane & Tetrahydrofuran		Isopropanol, n-Propanol		
2-Ethoxyethyl Acetate		Isopropyl Ether		
Acealdehyde		Isopropyl Myristate	-	
Acetic Acid 25%		Isopropylacetone		
Acetic Acid, Glacial		Kerosene		
Acetone		Magnesium Chloride aq. Sol.	-	
Ammonia, 25% ac. Sol.		Mercury		
Ammonium Hydroxide 25%		Methanol 98%		
Amyl & Propyl Acetate		Methyl Acetate		
Amyl Alcohol, Butanol		Methyl Cellosolve Acetate	-	
Aniline		Methyl Chloride	-	
Arsenic Acid		Methyl Ethyl Ketone	-	
Benzene		Methyle Isobutyl Ketone		
Benzyl Alcohol		Methylene Chloride		
Benzyl Benzoate		Monochlorobenzene, Freon		
Boric Acid 10%		Nitric Acid 25%	-	
Butyl Acetate		Nitric Acid 50%	-	
Calcium Chloride aq. Sol		Nitric Acid, Fuming		
Carbon Tetrachloride		Phenol, 100%		
Chlorine		Phenol, Aqueous 10%	-	
Chlorobenzene		Phosphoric Acid 25%		
Chloroform		Phosphoric Acid 85%		
Chromic Acid 20%		Potassium Chloride aq. Sol.		
Decalin		Potassium Hydroxide		
Diethylacetamide		Potassium Permanganate aq. Sol.		
Dimethyl Formamide		Propylene Glycol		
Dimethylsulphoxide (DMSO)		Propylene Glycol Acetate	_	
Ethanol 70%		Pyridine		
Ethanol 98%		Silicone Oil & Mineral Oil		
Ethyl Acetate		Silver Nitrate		
Ethyl Ether		Sodium Carbonate		
Ethylene Glycol		Sodium Dicromate		
Fluorinated Hydrocarbons		Sodium Hydroxide		
Formaldehyde Solution 30%		Sulphuric Acid 25%		
Formic Acid 25%		Sulphuric Acid, Concentrated		
Formic Acid 85%		Tetralin		
Gasoline		Toluene		
Glycerol		Trichloroacetic Acid 10%		
Hexane		Trichloroethylene		
Hydrochloric Acid 25%		Tricresyl Phosphate		
Hydrochloric Acid, Concentrated		Triethanolamine		
Hydrofluoric Acid 35%		Xylene		
Hydrogen Peroxide 30%		Zinc Chloride 10%		
lodine (tincture of)		Zinc Sulphate 10%		



6

Liquid Handling

Bottle Top Dispenser

Technical Information:

Product Code Volume Range	Volumo Pongo	Graduation	Accuracy		Coefficient of Variation	
	voluille hallye		%	μl	%	μl
550.001.305	0.5 – 5mL	0.1mL	0.5	25	0.1	5
550.001.310	1.0 – 10mL	0.2mL	0.5	50	0.1	10
550.001.325	2.5 – 25mL	0.5mL	0.5	125	0.1	25
550.001.350	5.0 – 50mL	1.0mL	0.5	250	0.1	50





Dispenser suitable for dispensing liquids with the following limits:

Temperature:	15 to 40 degrees Celsius	
Vapour Pressure:	Max. 500mBar	
Density:	Max. 2.2g/cm ³	

Ordering Information:

Each unit is supplied with: bottle top dispenser, 5 x PP adapters (GL25, GL28, GL32, GL38, S40), dismantling tool, filling tube (220mm) and user manual.

Product Code	Description	Pack Size
550.001.305	0.5 – 5mL	Each
550.001.310	1.0 – 10mL	Each
550.001.325	2.5 – 25mL	Each
550.001.350	5.0 – 50mL	Each

Accessories:

Product Code	Description	Pack Size
550.001.615	Tube Intake (PTFE) 15cm	Each
550.001.620	Tube Intake (PTFE) 20cm	Each
550.001.625	Tube Intake (PTFE) 25cm	Each
550.001.630	Tube Intake (PTFE) 30cm	Each



