

Vortex Shaker

Features & Benefits

- Small, touch or continuous tube mixer
- Vacuum suction feet hold the unit to the bench
- Aluminium cast base avoids unnecessary "walking"
- Versatile unit with a range of tube accessories

Technical Specifications:

Product Code	400.101.001
Voltage	230V/60Hz
Power	60W
Mixing Action	Orbital
Mixing Diameter	4mm
Motor Rating Input	58W
Motor Rating Output	10W
Speed Range	0 to 2,500rpm
Speed Display	Scale
Operation	Touch / Continuous
Dimensions	127 x 130 x 160mm
Weight	3.5kg
Permissible Ambient Temperature	5 to 40°C
Permissible Relative Humidity	80%
Protection Class	IP21



400.101.001 Vortex Shaker
Control Panel

Ordering Information:

Product Code	Description	Pack Size
400.101.001	Shaker Vortex supplied with 400.101.010 Standard Adapter	Each



400.101.001 Vortex Shaker
supplied with 400.101.010 Standard Adapter



400.101.001 Vortex Shaker
pictured with 50mL Centrifuge Tube
(Tube not supplied)

400.101.001 Vortex Shaker
pictured with 400.101.011 Rod &
400.101.015 Tube Adapter
(Micro tubes not supplied)

Ordering Information:

Product Code	Application	Pack Size
400.101.010	Standard Top for 30mm Tubes	Each
400.101.011	Tube Holding Rod for Tube Adapters	Each
400.101.012	Universal Adapter for Tube Adapters	Each
400.101.013	Tube Adapter for 6mm Ø x 48	Each
400.101.014	Tube Adapter for 10mm Ø x 18	Each
400.101.015	Tube Adapter for 12mm Ø x 12	Each
400.101.016	Tube Adapter for 16mm Ø x 8	Each
400.101.017	Tube Adapter for 20mm Ø x 8	Each
400.101.018	Platform Mat for Small Tubes & Vessels	Each
400.101.019	Rubber Suction Foot Replacement	Each



400.101.010
Standard Top
for 30mm Tubes



400.101.011
Tube Holding Rod for Tube Adapters



400.101.012
Universal Adapter
for Tube Adapters



400.101.013
Tube Adapter
for 6mm Ø x 48



400.101.014
Tube Adapter
for 10mm Ø x 18



400.101.015
Tube Adapter
for 12mm Ø x 12



400.101.016
Tube Adapter
for 16mm Ø x 8



400.101.017
Tube Adapter
for 20mm Ø x 8



400.101.018
Platform Mat for
Small Tubes & Vessels



400.101.019
Rubber Suction Foot
Replacement



400.101.015 Tube Adapter
fitted with 1.5mL Micro Tubes
(Micro tubes not supplied)



400.101.001 Vortex Shaker
fitted with 400.101.015 Tube adapter
& 400.101.012 Universal Attachment



400.101.001 Vortex Shaker
fitted with 400.101.012
Universal Attachment