Rotisserie Mixer

Features & Benefits

- Gentle but effective mix for samples up to 50mL
- Digital control for accuracy
- End over end design
- Wide variety of attachments available







400.104.040 with 400.104.050 (16 x 50mL) Attachment. Unit has ablility to be fully tilted and operated from 0 \leq 90°.



Technical Specifications:

| Product Code | 400.104.040 | |
|---------------------------------|--------------------|--|
| Voltage | 240V/60Hz | |
| Power | 40W | |
| Motor Type | DC Brushless | |
| Tilt Angle | 0 - 90° | |
| Speed Range | 10 to 70rpm | |
| Speed Display | Digital, LCD | |
| Timer | Yes | |
| Timer Display | Digital, LCD | |
| Timer Range | 1 to 1,199 Minutes | |
| Operation | Continuous / Timed | |
| Dimensions | 510 x 190 x 260mm | |
| Weight | 8kg | |
| Permissible Ambient Temperature | 5 to 40°C | |
| Permissible Relative Humidity | 80% | |
| Protection Class | IP21 | |



400.104.040 Display/Control Panel

Ordering Information:

| Product Code | Description | Pack Size |
|--------------|---|-----------|
| 400.104.040 | Rotisserie Mixer supplied with 400.104.050 Platform | Each |

Rotisserie Mixer - Accessories

Ordering Information:

| Product Code | Application | Pack Size |
|--------------|-----------------------------|-----------|
| 400.104.045 | For 1.5mL x 48 (Horizontal) | Each |
| 400.104.046 | For 15mL x 24 (Horizontal) | Each |
| 400.104.047 | For 50mL x 24 (Horizontal) | Each |
| 400.104.048 | For 1.5mL x 32 (Vertical) | Each |
| 400.104.049 | For 15mL x 16 (Vertical) | Each |
| 400.104.050 | For 50mL x 16 (Vertical) | Each |





400.104.045 Horizontal Attachment for 1.5mL x 48 Micro Tubes (tubes not included)



400.104.048 Vertical Attachment for 1.5mL x 32 Micro Tubes (tubes not included)



400.104.046 Horizontal Attachment for 15mL x 24 Centrifuge Tubes (tubes not included)



400.104.049 Vertical Attachment for 15mL x 16 Centrifuge Tubes (tubes not included)



400.104.047 Horizontal Attachment for 50mL x 24 Centrifuge Tubes (tubes not included)



400.104.050 Vertical Attachment for 50mL x 16 Centrifuge Tubes (tubes not included)